

Required Records for Surgical Guides

Tissue Supported

Bone Supported

With Provisional

Without Provisional

With Provisional

Without Provisional

Fully Dentate

- CBCT Scan of Patient
- Upper and Lower Impressions
- Bite Registration
- Photos

Fully Dentate

- CBCT Scan of Patient
- Upper and Lower Impressions
- Bite Registration

Fully Dentate

- CBCT Scan of Patient
- Upper and Lower Impressions
- Bite Registration
- Photos

Fully Dentate

- CBCT Scan of Patient
- Working Arch Impression

Partially Dentate

- CBCT Scan of Patient Wearing Wax Rim with Markers
- Send Physical Wax Rims
- Upper and Lower Impressions
- Bite Registration
- Photos

Partially Dentate

- CBCT Scan of Patient Wearing Wax Rim with Markers
- Send Physical Wax Rims
- Upper and Lower Impressions
- Bite Registration

Partially Dentate

- CBCT Scan of Patient Wearing Wax Rim with Markers
- Send Physical Wax Rims
- Upper and Lower Impressions
- Bite Registration
- Photos

Partially Dentate

- CBCT Scan of Patient of Patient Wearing Wax Rim with Markers
- Send Physical Wax Rims
- Working Arch Impression

Fully Edentulous

- CBCT Scan of Patient Wearing Denture or Wax Rim
- Send Physical Wax Rims (if applicable)
- Upper and Lower Impressions
- Bite Registration
- Photos

Fully Edentulous

- CBCT Scan of Patient Wearing Denture or Wax Rim
- Send Physical Wax Rims (if applicable)
- Upper and Lower Impressions
- Bite Registration

Fully Edentulous

- CBCT Scan of Patient Wearing Denture or Wax Rim
- Send Physical Wax Rims (if applicable)
- Upper and Lower Impressions
- Bite Registration
- Photos

Fully Edentulous

- CBCT Scan of Patient Wearing Denture or Wax Rim
- Send Physical Wax Rims (if applicable)
- Working Arch Impression

Required Records for Surgical Guides

Teeth and Tissue

Teeth and Bone

With Provisional

Without Provisional

With Provisional

Without Provisional

Fully Dentate

- CBCT Scan of Patient
- Upper and Lower Impressions
- Bite registration
- Photos

Fully Dentate

- CBCT Scan of Patient
- Upper and Lower Impressions
- Bite registration

Fully Dentate

- CBCT Scan of Patient
- Upper and Lower Impressions
- Bite registration
- Photos

Fully Dentate

- CBCT Scan of Patient - Working Arch Impression

Partially Dentate

- CBCT Scan of Patient Wearing Wax Rim with Markers
- Send Wax Rim with Markers
- Upper and Lower Impressions
- Bite registration
- Photos

Partially Dentate

- CBCT Scan of Patient Wearing Wax Rim with Markers
- Send Wax Rim with Markers
- Upper and Lower Impressions
- Bite registration

Partially Dentate

- CBCT Scan of Patient Wearing Wax Rim with Markers
- Send Wax Rim with Markers
- Upper and Lower Impressions
- Bite Registration
- Photos

Partially Dentate

- CBCT Scan of Patient Wearing Wax Rim with Markers
- Send Wax Rim with Markers
- Working Arch Impressions